

**We bring quality
to light.**



Measurement of Near Eye Display Using Goniophotometer and Imaging Light Measuring Device

Tianxing Zhu | 2025

Your Application – our Know-how!

Display Measurement



AR/VR



Automotive
Interior Lighting



Display
Characterization



Flat Panel Display
Production




MicroLEDs



System Audit

Light Emitter Measurement



LED Production



LED / Modules



UV-LEDs



VCSEL / Laser



IR Sensing

Lighting Products Measurement



Automotive
Exterior Lighting



Blue Light
Hazard



Solid State
Lighting

AR/VR Testing Components



Final Headset



Alignment



Display-Lens Module
& Optics



Microdisplay




Final Headset




Prescription



Waveguide



Microdisplay / Projector



Agenda

- 1. Introduction**
- 2. AR/VR Near-Eye-Display Measurement Systems**
- 3. Measurement Setup**
- 4. Result Comparison**
- 5. Conclusions**

Introduction

The quality/characteristic of a near eye display (NED) is important for the final performance of the VR/AR headset.

In this presentation two different measurement setups for near eye displays are compared:

- ▲ Goniophotometer DMS 803 (Instrument Systems)
- ▲ Imaging light measurement device (ILMD) LumiTop 5300 ARVR (Instrument Systems).

AR/VR Near-Eye-Display Measurement



DMS Goniophotometer

- ✓ Angle resolved spectral spot measurement
- ✓ Eyebox / Qualified Viewing space
- ✓ Time resolved




CAS spectrometer as reference

- Highest Spectral resolution
- Absolute calibration of luminance and color
- Traceable accuracy $Y = \pm 3\%$, $x/y = 0.0015$



TOP 300 AR/VR Optical Probe

- ✓ Robust, compact and lightweight design for manufacturing applications
- ✓ High-accuracy luminance and color measurement
- ✓ Integrated viewfinder camera with 5 MP
- ✓ Optics mimic the human eye with a FOV of $\pm 1.2^\circ$
- ✓ Binocular setup possible





LumiTop AR/VR Imaging Light Measurement Device

- ✓ 2D one-shot measurement at 120° FOV
- ✓ Distortion, Sharpness, Uniformity, Contrast, ...
- ✓ Time resolved
- ✓ LumiSuite for integration and analysis

DMS Goniophotometer – Reference System

DMS-803: Reference Goniophotometer


1. 6 motorized axes for automated measurements
 - Fully automated spectral sampling
 - Dedicated optics to mimic human visual system captures both: eye rotation and retinal image
 - Variety of accessories for ambient, specular and backlights are suitable for your angular measurements.
2. CAS 140D reference spectroradiometer
 - Spectral spot information for accurate measurement values
 - High-end CCD detector cooled and back illuminated




LumiTop AR/VR

LumiTop: 3-way beamsplitter (→ one-shot)


1. Camera system for 2D measurement
 - Calibrated System
 - 120° FOV & 40 px/deg.
2. CAS 140D reference spectroradiometer
 - Spectral spot information for live color fine-correction
 - High-end CCD detector cooled and back illuminated
3. Flicker diode
 - Sync to display modulation
 - Enables High Luminance mode



Differences between Spot and 2D Measurements



Luminance as seen by DMS angular resolved goniometric spot measurements.



Luminance as seen by LumiTOP spectrally enhanced imaging colorimeter.


Measurements

▲ Several characteristics are measured to obtain a perfect VR device


- Final device
 - Luminance
 - Chromaticity
 - Gamma
 - Contrast
- Display
 - Defects
 - Mura
 - Luminance
 - Chromaticity
 - Gamma
 - Contrast
- Optical module (display+optics)
 - Distortion
 - MTF
 - FOV
 - VID
 - Chromaticity
 - Black/White
 - Grey 255
 - Contrast
 - Uniformity (Luminance/Color)
 - Chromatic aberration
 - Eye box
 - Veiling glare index

Measurement Setup

▲ DMS-803 Setup





▲ LumiTop 5300 AR/VR



DUT: Pancake + LCD; 90° FOV


Alignment

- ▲ Proper alignment is needed for comparison.
- ▲ E.g. a test pattern with crosses or blobs are used to check orientation and positioning.







Color comparison DMS/LT

- ▲ The sample was measured with DMS and with LumiTop and the measurements are compared.
- ▲ The pixels of the LumiTop are averaged to fit the FOV of the DMS.
- ▲ DMS settings
 - Azimuth 0°-355°; Step size 5°
 - Elevation 0°-50°; Step size 2°
 - 1872 Measurement points
- ▲ LumiTop covers complete FOV **with one image (no need to move)**
- ▲ For illustration: Camera image overlayed with the measurement points of the DMS (black circles)



Color comparison DMS/LT - white



- DMS-white: $\Delta u'v'$ relative to center
- LT-white: $\Delta u'v'$ relative to center
- DMS-LT-white: $\Delta u'v'$ between DMS-LT
- Good match between DMS and LT measurements



Comparison LumiTop AR/VR :: DMS ARVR

- Display measurement with Pancake optics
 - DUT NED with LCD and pancake lens

Difference between LMD and ILMD (DMS-LT) for different color pattern	Mean $\Delta u'v'$	Standard deviation $\Delta u'v'$
White	0.0005	0.0006
Red	0.0016	0.0010
Green	0.0028	0.0025
Blue	0.0016	0.0012



Conclusions

- ▲ The characterization of NED is done with an LMD and an ILMD.
 - Absolute, accurate and traceable calibration of light measurement devices enables comparability between different measurement principles.
 - No DUT specific calibration (“golden sample”) is required, even across different display technology platforms, with absolute calibration of the measurement device.
- ▲ We were able to quantitatively identify correlations between one-shot wide FOV and goniometric narrow FOV MTF measurement.
 - These correlations, based on high resolution goniometric reference measurements, could be used to speed up quality control measurements using wide FOV one-shot MTF measurements.

THANK YOU

for your attention!